


+ OneNPWT clinical decision tree for open wounds



The PICO pumps contain a MAGNET. Keep the PICO pumps at least 4 inches (10cm) away from other medical devices at all times. As with all electrical medical equipment, failure to maintain appropriate distance may disrupt the operation of nearby medical devices[†]. For detailed product information, including indications for use, contraindications, effects, precautions, warnings, and important safety information, please always consult product's user manual instructions for use (IFU) prior to use.

The information herein is intended for healthcare professionals. RENASYS is contraindicated in the presence of untreated osteomyelitis, exposed arteries/veins/organs/nerves, necrotic tissue with eschar present, malignancy in the wound, non-enteric and unexplored fistulas, and exposed anastomotic sites. Excessive bleeding is a serious risk associated with the application of suction to wounds, which may result in death or serious injury. For full product and safety information, please see the Instructions for Use.

* Wounds must not contain exposed arteries, veins, nerves or organs † Wounds should generally be no more than 4.5cm in depth and must not contain exposed arteries, veins, nerves or organs ‡ p=0.046; n=31; Compared to black foam in acute post traumatic wounds. **Reference: 1.** Brownhill R. PICO[®] Biomechanical Study. Data on file report. August 2019. DS/19/211/R. **2.** Hurd T, Chadwick P, Cote J, Cockwill J, Mole T, Smith J. Impact of gauze-based NPWT on the patient and nursing experience in the treatment of challenging wounds. International Wound Journal. 2010;7(6):448-455. **3.** Fraccalvieri M, Scalise A, Ruka E, et al. Negative pressure wound therapy using gauze and foam: Histological, immunohistochemical, and ultrasonography morphological analysis of granulation and scar tissues - Second phase of a clinical study. In: European Journal of Plastic Surgery. Vol 37 2014:411-416. **4.** Johnson S. V1STA[®] - A new option in Negative Pressure Therapy. Journal of Wound Technology. 2008;1:30-31. **5.** Fraccalvieri M, Ruka E, Bocchiotti M, Zingarelli E, Bruschi S. Patient's pain feedback using negative pressure wound therapy with foam and gauze. International wound journal. 2011;8(5):492-499. **6.** Smith+Nephew 2009. A prospective, open labelled, multicentre evaluation of the use of VISTA in the management of chronic and surgical wounds and A prospective, open labelled evaluation of the use of EZCare in the management of chronic and acute wounds. Internal Report. SR/CIME/010/012. [◇]Trademark of Smith+Nephew. All Trademarks acknowledged. ©2021 Smith+Nephew. 28194_US